



Viral RNA Extraction Kit (Spin Column)

Intended Use

Used for nucleic acid extraction with high nucleic acid purity in subsequent experiments (DNA/RNA)

Product Advantages



The whole extraction time is less than 20min, simple and fast, suitable for nucleic acid extraction of most biological samples.



The product can be transported at room temperature to maintain enzyme activity and reduce transportation cost.



High purity nucleic acid can be obtained as well as high yield by using silica gel membrane as substrate.



Ordering Information

Catalog #	Kit Size	Suitable Sample	Main Components
P082T00401	24 Tests/Kit	Stool, Swab, Sputum, Serum and Plasma	Lysis buffer Washing buffer A&B
P082T00501	50 Tests/Kit		Elution buffer Spin Column
P082T00601	100 Tests/Kit		Centrifugal tubes Carrier RNA * Proteinase K *

* Carrier and Proteinase K are optional for Stool, Swab and Sputum sample, please contact our sales for more detailed information.



Jiangsu Mole Bioscience Co.,Ltd

6-7 th Floor, C116 Building, No.805, Jiankang Avenue,
Medical New&Hi-tech District, Taizhou,
Jiangsu Province, P.R.C.

Web: www.molechina.com
Email: info@molechina.com
Tel: +86-571- 87209310