


6 Minutes

High Purity

Double Protection

NEW UPGRADE

6 Minutes MOLE Fast 96 Nucleic Acid Extraction Instrument



Safe

Accurate

Fast

Product Features

● Fast Extraction

Fast extraction in 6 minutes.

● Accurate Results

High purity extraction, magnetic beads recovery rate $\geq 99.5\%$.

● Double Protection

Ultraviolet sterilization device and 11-level high-efficiency air filtration system can effectively reduce the aerosol pollution.

● Power-off Protection

Program keep working after power interruption.

● Liquid Volume Identification

The amplitude can be automatically adjusted according to the liquid volume, avoid liquid overflow due to improper parameter settings.

● Processing Capacity

The program is compatible with 96 and 24 samples extraction, and whole process will be finished in 6 minutes.

Product Parameter

Technical Parameters	96 Wells	24 Wells
Processing Capacity	96	24
Processing Volume	20-1000 μ l	100 μ l-5ml
Model No.	Mole Fast96	
Magnetic Beads Recovery Rate	$\geq 99.5\%$	
Power	AC100-240V, 50-60Hz	
Screen	7 Inches, 870*1550 Color Screen in Chinese or English	
Weight	45kg \pm 0.5kg	
Size	(685mm*840*390) mm \pm 10 mm	
UV Lamp	Equipped	
Pollution Control	Closed cover and 11-Level High-Efficiency Air Filtration System	



Jiangsu Mole Bioscience Co., Ltd.
Stock code: 839185
info@molechina.com

Jiangsu-Taizhou
Add: 6-7th Floor, G116 Building, No.805, Jiankang Avenue, Medical New&Hi-tech District, Taizhou, Jiangsu Province, China
Tel: 0523-86080606

Zhejiang-Hangzhou
Add: 4-5th Floor, Building No.5, Health Valley, No.1500 Wenyi xi Road, Yuhang District, Hangzhou, Zhejiang province, China
Tel: 0571-87209306 / 87209308