



Viral DNA Extraction Kit (Boiling Method)

Intended Use

This kit is intended for isolation of DNA of viruses, bacteria or human cells from feces, swabs, gastric mucosa and cervical cytological specimens.



Product Advantages



The product can be used for rapid nucleic acid extraction in large quantities (DNA)



The product can be transported at room temperature to maintain enzyme activity and reduce transportation cost.



The product does not involve toxic organic solvents, such as phenol, chloroform, etc.

Ordering Information

Catalog #	Kit Size	Suitable Sample	Main Components
P031T020A0C0	20 Tests/Kit	Feces, swabs, gastric mucosa tissue, cervical cytological specimens.	Lysis buffer (EDTA, NP-40, Chelex-100)
P031T024A0C0	24 Tests/Kit		
P031T050A0C0	50 Tests/Kit		
P031T100A0C0	100 Tests/Kit		

Jiangsu Mole Bioscience Co., Ltd.
6-7th Floor, G116 Building, No.805,
Jiankang Avenue, Medical New&Hi-tech District,
Taizhou, Jiangsu Province, China

Tel: 0571-87209310
Email: info@molechina.com
www.molechina.com