


CE IVD




# COVID-19 and Influenza A & B Rapid Test Cassette


Qualitative membrane strip-based immunoassay for the detection of influenza A virus, influenza B virus, and COVID-19 virus.




## Kit Components




Test Cassette



Extraction buffer



Swab



Workstation



Instruction for use

Jiangsu Mole Bioscience Co., Ltd.  
6-7th Floor, G116 Building, No.805,  
Jiankang Avenue, Medical New&Hi-tech District,  
Taizhou, Jiangsu Province, China

Tel: 0571-87209310  
Email: info@molechina.com  
www.molechina.com

## Product Specification

### SAMPLE TYPES

Nasopharyngeal Swab  
Nasal Swab

### SHELF TIME

24 Months

### STORAGE CONDITION

2-30 °C

## Product Features



Cost Effective



User Friendly



Fast Result in 15 mins


#### All in 1 Testing

COVID-19, Influenza A & B can be detected and differentiated by one sampling.

#### Better Performance

The sensitivity and specificity for Influenza A or B are both >95%.  
The sensitivity and specificity for COVID-19 are both >98%.

## Test Procedure



## Ordering Information

Cat No.	Product Description	Kit Size	Certificate
R218T020B0C0	Flu A/B + COVID-19 Antigen Rapid Test ( 2 Windows)	20 Tests/Kit	CE
R218T020B1C0	Flu A/B + COVID-19 Antigen Rapid Test ( 3 in 1 Window)	20 Tests/Kit	CE
R501T020B0C0	Flu A+ B Antigen Rapid Test ( 2 in 1 Window)	20 Tests/Kit	CE