

PURE 32



Automatic Nucleic Acid Extractor



Cat No.: MLN32P201



Plug Type: European Standard Plug



Sample numbers.: 1-32



LCD: 480x800 Color display



Weight: 23±1 kg



Processing volume: 20-1000μL



Temperature control: 10-99°C



Power: AC 110-240 V, 50-60 Hz



Magnetic bead recovery rate: ≥98%



Product Size: 350x380x415 mm



Product features



Advanced Oscillation Technology



Safety & Reliability



Contamination Control



High Magnetic Bead Recovery



Jiangsu Mole Bioscience Co., Ltd.
6-7th Floor, G116 Building, No.805,
Jiankang Avenue, Medical New&Hi-tech District,
Taizhou, Jiangsu Province, China

Tel: 0571-87209310
Email: info@molechina.com
www.molechina.com

		Applicable specimen types										Extraction time	Kit size
	Swabs	Others											
Cat No.	Nasopharynx, Oropharynx, Nasal and Anal	Sputum	Stool	Urine	Whole Blood	Serum/ Plasma	Herpes Exudate	Cervical Shedding Cell	FFPE				
P083/P084	●	●	●							13 mins			
P085					●					30 mins			
P086						●				24 mins			
P088				●		●	●	●		20 mins			
P089			●							31 mins			
P090	●	●	●	●	●	●	●	●		31 mins			
P091	●	●	●							31 mins			
P120	●	●				●				30 mins			
P140									●	21 mins			