

Clarithromycin, Quinolone-Resistant and *Helicobacter Pylori* Nucleic Acid Amplification Test Kit

Fluorescent Probe-Based Real-Time PCR Assay

Why do we test *Helicobacter Pylori*?

About 50% of people in the world have a history of *Helicobacter Pylori* infection; The Maastricht V/Florence Consensus Report pointed out “*H. pylori* infection is accepted as the major aetiological factor for gastric cancer, and *H. pylori* eradication reduces the risk of gastric cancer development.”

Why is it important to detect *Helicobacter Pylori* antibiotic resistance?

Due to the extensive use of antibiotics, *H. pylori* resistance strains are widespread; The Maastricht V/Florence Consensus Report recommended performing clarithromycin susceptibility testing by the molecular test directly.

Test Principle

This kit uses polymerase chain reaction (PCR) based on TaqMan fluorescent probe technology to identify *H. pylori* and antibiotic-resistant mutations.



Clinical significance

- Select** effective antibiotics and treatment strategies to eradicate *H. pylori*.
- Reduce** the cost of diagnosis and treatment of *H. pylori* infection.
- Avoid** dysbiosis due to multiple antibiotics and prolonged treatment duration.
- Prevent** future antibiotic resistance due to the failure of first line treatment.


Comparison of eradication rate with susceptible strains and resistant strains after 14 days of treatment.

Treatment protocol	Susceptible strains	CAM and/or LVX Resistant strains
Standard triple therapy Levofloxacin triple therapy Bismuth quadruple therapy	Eradication rate > 95%	Eradication rate < 50%

Recommended Thermal Cyclers

- Applied Biosystems™ 7500 Real-Time PCR System
- SLAN 96 Real-Time PCR System
- Other Real-Time Thermal Cyclers which could detect FAM, VIC/HEX, ROX and CY5 simultaneously

Sample Results



Ordering Information

Cat No.	Product Discription	Format	Specimen	Shelf Life	Kit Size	Certification
B052T050B0C2	Mospire Clarithromycin, Quinolone-resistant and Non -resistant <i>Helicobacter Pylori</i> Nucleic Acid Amplification Test Kit (Fluorescent Probe-Based Real-Time PCR Assay)	Liquid	Gastric biopsy	12 Months	50 Tests/Box	CE NMPA

F-B052-001-EN Ver: 2.1 Date: 2026-01-22



Jiangsu Mole Bioscience Co., Ltd.
No.1, Lingzhi Road, Medical New & Hi-tech District (Gaogang District), 225300 Taizhou, Jiangsu Province, China

Tel: 0571-87209310
Email: info@molechina.com
www.molechina.com