

SARS-COV-2

ANTIBODY (IgG & IgM) RAPID TESTS

- Antibodies are proteins created by your body's immune system soon after you have been infected or vaccinated.
- Easy to use, no instrument and location required.
- Rapid Result in 15 min.
- Detection of different stages of virus infection.

Why Using Antibody Rapid Test?



Test Cassette



Dropper



Buffer




Package Insert



KIT COMPONENT

Ref	Product Name	Packing Size	Storage Condition	Certification
R206T020B0	A rapid IgM-IgG combined antibody Test kit for SARS-CoV-2	20 Tests/Kit	2-30 °C	CE

RESULTS INTERPRETATION



Jiangsu Mole Bioscience Co., Ltd.
Stock code: 839185
info@molechina.com

Jiangsu-Taizhou
Add: 6-7th Floor, G116 Building, No.805, Jiankang Avenue, Medical New&Hi-tech District, Taizhou, Jiangsu Province, China
Tel: 0523-86080606

Zhejiang-Hangzhou
Add: 4-5th Floor, Building No.5, Health Valley, No.1500 Wenxi Road, Yuhang District, Hangzhou, Zhejiang province, China
Tel: 0571-87209306 / 87209308