

Saliva Collection Kit

Simple • Flexible • Convenient • Safe • Extensive



The Saliva specimen Collection Kit is suitable for collecting and storing saliva specimen at room temperature. It is widely used in DNA gene detection, cytological sample processing, virology cell culture, enzyme immunoassay (ELISA) detection, PCR detection, molecular-based diagnosis detection and saliva extraction of microorganisms.

Product Advantages

- **Simple:** The collection process is simple and convenient, and can be operated by individual.
- **Flexible:** The sampling site is not limited by space and site, and it can realize contactless sampling at home.
- **Convenient:** Saliva specimen are stable, can be stored at room temperature and easy to transport.
- **Safe:** Noninvasive sampling reduces the risk of infection among health care workers.
- **Extensive:** Suitable for all kinds of people, achieving rapid collection of specimens.

Product Parameters

Packing specification:	1 set/kit, 10 sets/kit	Kit content:	5 mL Vial with a Funnel, a Separate 5 mL Vial Filled with 2 mL buffer, label sticker and biosafety bag
Storage condition:	Room temperature (2-30°C)	Expiration:	36 months

江苏默乐生物科技股份有限公司
Jiangsu Mole Bioscience Co., Ltd.

Stock code : 839185

Jiangsu-Taizhou
Add: 6-7th Floor, G116 Building, No.805, Jiankang Avenue, Medical New&Hi-tech District, Taizhou, Jiangsu Province, China
Tel: 0523-86080606

Zhejiang-Hangzhou
Add: 4-5th Floor, Building No.5, Health Valley, No.1500 Wenyixi Road, Yuhang District, Hangzhou, Zhejiang province, China
Tel: 0571-87209306 / 87209308