

PURE 32



Automatic Nucleic Acid Extractor



Cat No.: Pure32



Plug Type: European Standard Plug



Sample numbers.: 1-32



LCD: 480x800 Color display



Weight: 23±1 kg



Processing volume: 20-1000μL



Temperature control: 10-99°C



Power: AC 110-240 V, 50-60 Hz



Magnetic bead recovery rate: ≥98%



Product Size: 350x380x415 mm



Product features



Advanced Oscillation Technology



Safety & Reliability



Contamination Control



High Magnetic Bead Recovery



Jiangsu Mole Bioscience Co., Ltd.
6-7th Floor, G116 Building, No.805,
Jiankang Avenue, Medical New&Hi-tech District,
Taizhou, Jiangsu Province, China

Tel: 0571-87209310
Email: info@molechina.com
www.molechina.com

	Applicable specimen types									Extraction time	Elution buffer volume	Exp. of downstream testing	kit size
	Swabs	Others											
Cat No.	Nasopharynx, Oropharynx, Nasal and Anal	Sputum	Stool	Urine	Whole Blood	Serum/ Plasma	Herpes Exudate	Cervical Shedding Cell	FFPE				
P083/P084	●	●	●							13 mins	100μL	IFA, IFB	16 tests/plate, 1 test/strip (Single package).
P085					●					30 mins	100μL	EB, B27	
P086						●				24 mins	50μL	HBV, Dengue	
P088				●		●	●	●		20 mins	50μL	HPV	
PL088				●		●	●	●		20 mins	200μL	HPV	
P089			●							31 mins	500μL	HP (Stool)	
P090	●	●	●	●	●	●	●	●		31 mins	400μL	META Urogenital Panel	
P091	●	●	●							31 mins	400/230μL	META Respiratory Panel	
P120	●	●				●				30 mins	50μL	TB	
P140									●	21 mins	50μL	K-Ras, HP (FFPE Gastric Biopsy)	



