











# PURE 32



## Automatic Nucleic Acid Extractor

-  Cat No.: Pure32
-  Plug Type: European Standard Plug
-  Sample numbers.: 1-32
-  LCD: 480x800 Color display
-  Weight: 23±1 kg
-  Processing volume: 20-1000μL
-  Temperature control: 10-99°C
-  Power: AC 110-240 V, 50-60 Hz
-  Magnetic bead recovery rate: ≥98%
-  Product Size: 350x380x415 mm



### Product features



Advanced Oscillation  
Technology



Safety & Reliability



Contamination Control



High Magnetic Bead  
Recovery


Ver: 3.0 Date: 2024-05-06



Jiangsu Mole Bioscience Co., Ltd.  
No.1, Lingzhi Road, Medical New & Hi-tech  
District (Gaogang District), 225300 Taizhou,  
Jiangsu Province, China

Tel: 0571-87209310  
Email: info@molechina.com  
www.molechina.com

## Ordering Information



### Applicable specimen types

Cat No.	Swabs		Others								Extraction time	Elution buffer volume	Exp. of downstream testing
	Nasopharynx, Oropharynx, Nasal and Anal	BAL* NAL*	Sputum	Stool	Urine	Whole Blood	Serum/ Plasma	Herpes Exudate	Cervical Shedding Cell	FFPE			
P083/P084	●		●	●							13 mins	100μL	IFA, IFB
P085						●					30 mins	100μL	EB, B27
P086							●				24 mins	50μL	HBV, Dengue
P088					●		●	●	●		20 mins	50μL	HPV 2+16
PL088					●		●	●	●		20 mins	200μL	HPV 24 genotyping
P089				●							31 mins	500μL	HP (Stool)
P090	●	●	●		●	●	●	●	●		31 mins	400μL	META Urogenital Panel
P091	●	●	●								31 mins	400/230μL	META Respiratory Panel
P120	●		●				●				30 mins	50μL	TB
P140										●	21 mins	50μL	K-Ras, HP (FFPE Gastric Biopsy)

\*BAL: Bronchoalveolar lavage fluid

\*NAL: Nasal lavage fluid

Kit Size: Both Plate (16 tests/plate) and Single Package (1 test/strip) are available.

info@molechina.com | www.molechina.com